

OPENLAB VERIFIED · Report identity and source-transcribed quantity. Compound purity and analytical events retain their separate source scopes.

M-STEN

Lab record

Batch	Laboratory	Report	Analysed
OL2208	Janoshik Analytical	28941	18 April 2023

Powered by EvidenceOS, connecting this Olympus Labs UK product and batch to its laboratory results and source evidence.

M-STEN · Concentration

7.05 mg

Reported sample; the laboratory gives mg without a per-capsule denominator.

Janoshik report 28941 · OL2208

M-STEN · Purity

99%

Compound raw material

Olympus Labs UK · compound purity statement

REPORTED

Compound identity

Methylstenbolone

REPORTED

Reported quantity

7.05 mg

REPORTED

Compound purity

99%

REPORTED

Laboratory source

Janoshik

REPORTED

Batch connection

OL2208

REPORTED

Sample receipt

2023-04-14

The batch report and supporting compound analyses retain their separate dates and source identities. Compound raw-material purity is a supporting observation; production stage is unspecified.



Compound analysis

Methylstenbolone · LC/MS chromatogram · 2022-08-16

LC/MS 16 August 2022 DAD1 B, 220 nm

LC/MS Report

File D:\CHEM32\DATA\202208\20220816-KEHU-1 2022-08-16 13-53-55\MS-20220816-LCMS.D
Injection Date : 16-Aug-22, 14:22:09 Tgt Mass (EZ) :
Sample Name : Location : P1-F-08
Acq. Operator : Inj : 1
Spec. Reported : MS Integration Inj Volume : 6 ul
Acq. Method : D:\CHEM32\DATA\202208\20220816-KEHU-1 2022-08-16 13-53-55\401-CZ-RK.M
Analysis Method : C:\Chem32\1\METHODS\401-GP50-KEHU.M
Sample Info :
Method Info : Instrument : Agilent 1200-6110 MS
Mobile Phase : A=0.1%TFA in water, B=Acetonitrile
Gradient : 5-95% B in 4 min, 95%B in isocratic for 3min
Flow Rate : 0.6 mL/min
Column : XDB-C18 2.1*50mm, 3.5um
Column temperature : 35°C



Integration Results for DAD1 B, Sig=220,8 Ref=off

RetTim	Width	Area	Height	Area%
13.18	0.04	9.47	3.49	0.09
13.23	0.05	15.30	4.21	0.15
13.30	0.08	23.18	4.86	0.22
14.16	0.05	9.13	3.29	0.09
14.26	0.11	30.66	4.80	0.29
14.84	0.08	3378.06	673.75	32.43
14.94	0.10	6914.67	1016.36	66.38
15.21	0.07	16.34	3.20	0.16
15.31	0.08	19.68	3.78	0.19


Compound analysis

Methylstenbolone · mass spectrum · 2022-08-16


LC/MS 16 August 2022 Original source spectrum

LC/MS Report

Ret. Time: 15.00 <<<< POSITIVE SPECTRA >>>>



Ret. Time: 15.43



Original analytical page · Supporting compound evidence

Sample analysis

M-STEN · PDA Ch1 220–260 nm · 2023-04-14

Chromatography 14 April 2023 PDA Ch1 220–260 nm




Analysis Report

<Sample Information>

Sample Name : 1275-4
Method Filename : Raws - blind common anabolic steroid and aromatase inhibitors screening.lcm
Vial # : 3-4
Injection Volume : 0.5 uL
Date Acquired : 14/04/2023 10:41:00
Date Processed : 14/04/2023 10:45:35
Sample Type : Unknown
Acquired by : Janoshik s.r.o.
Processed by : Janoshik s.r.o.

<Chromatogram>

mAU



<Peak Table>

PDA Ch1 220-260nm

Peak#	Ret. Time	Area	Name
1	1.600	114362	
Total			

info@janoshik.com

Original analytical page · Report 28941

Sources and record edition

TEST REPORT



E-mail: info@janoshik.com
Web: www.janoshik.com

Task Number

#28941

Testing ordered >

31 MAR '23

Sample received >

14 APR '23

Client

support@olympuslabs.uk

Sample

M-STEN 6.6MG

Manufacturer

Olympus Labs UK

Batch

OL2208

Sample description >

Pictures not included.

Tests requested >

Qualitative and quantitative analysis of an AAS sample.

Results >

Methylstenbolone	7.05 mg

Comments >

Analysis conducted >

18 APR 2023

Signature >



Verify this test at www.janoshik.com/verify/ with the following unique key

Z5TLN9FNRFJ2

Laboratory dates

Sample received: 14 April 2023

Analysed: 18 April 2023

Janoshik · 28941

[View source](#)

25444396ece3afd08e492391ae5817e33b5e1a1c902a5c5a558001c1f3bfdcce

Methylstenbolone · LC/MS chromatogram · 2022-08-16

Source image included in this edition.

2d5d05b546b92962528ef10fa9f2bf08167a4dc4c4a0ce9e2dd865ed925e5810

Methylstenbolone · mass spectrum · 2022-08-16

Source image included in this edition.

ce6973dfa3086df54711e5b2f3c95e98dfaaf9ce98009c01d1215f049ef6b1f1

Janoshik report 28941

[View source](#)

ab3b9e45833e5fddb213bf048e9575351930bff00ec3a1ef22f7c00158dd9099

M-STEN · publisher purity statement

[View source](#)

b5492c937ef156a81746323e2fc7f95b48611c9b145990b1eec473ad4110b42b

M-STEN · PDA Ch1 220-260 nm · 2023-04-14

[View source](#)

cbf09787f1fcefef2503a3798b4131c9c6eed3ab208dfc145460e4fc55a517212

Edition

openlab.v2.9df4b788bb0bc0e4

Record digest:

8ff76b0755ad06084e50cf49aa919349c3be61ca8d5d60a63291f246927ad98a